

Catalog Number

75020-SP50

Semi Prep HPLC Column, HPS Phenyl, 5 μ m, 120A. 10mm ID x 50mm Length. Cogent Brand. 1 EA.



Additional Info

This Cogent HPS Phenyl™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is fully end capped and is manufactured with a mono functional, 4 carbon ligand that has an attached *Phenyl* group and it is bonded to the highest purity, 120A silica available.

With a 5 μ m Particle, this Column works well with higher flow rates producing lower back pressures compared to columns with smaller particles making ideal for the purification of Isomers (*positional*), Tocopherols, Flavonoids, Polynuclear Aromatics, Nitro-Aromatic Compounds, Active Pharmaceutical Ingredients (APIs), Natural

Price

\$2,588.73

Additional Info

Products and Related Compounds in Reversed Phase HPLC.

Compare this column's performance to Luna™ or Inertsil™ columns.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$2,588.73