

# Catalog Number

# 75025-P250

**Prep HPLC Column, HPS Cyano, 5 $\mu$ m, 120A. 21.2mm ID x 250mm Length. Cogent Brand. 1 EA.**



## Additional Info

This Cogent HPS Cyano™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is fully end capped and is manufactured with a mono functional, 3 carbon ligand that has an attached terminal Cyano (*Carbon & Nitrogen*) group which includes a triple bond, and it is bonded to our 120A High Purity Silica, making this Column ideal for purification of polar compounds in Reversed Phase.

With a 5  $\mu$ m Particle, this Column works well with higher flow rates and lower back pressures while the large diameter and length of this column provides high loading

**Price**

**\$13,119.21**

## Additional Info

capacity which is great for Preparative HPLC.

*Compare the performance of this Column to Inertsil or Luna brands.*

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$13,119.21**