

Catalog Number

75025-SP250

Semi Prep HPLC Column, HPS Cyano, 5 μ m, 120A. 10mm ID x 250mm Length. Cogent Brand. 1 EA.



Additional Info

This Cogent HPS Cyano™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is fully end capped and is manufactured with a mono functional, 3 carbon ligand that has an attached terminal Cyano (*Carbon & Nitrogen*) group which includes a triple bond, and it is bonded to our 120A High Purity Silica, making this Column ideal for purification of polar compounds in Reversed Phase.

With a 5 μ m Particle, this Column works well with higher flow rates and lower back pressures while the large diameter and length of this column provides high loading

Additional Info

capacity which is great for Preparative HPLC.

Compare the performance of this Column to Inertsil or Luna brands.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com