

Catalog Number

82000-SP-GH

Guard Column Holder, Semi Prep, 10mm TYPE-C Guard. Accepts 10mm ID x 10mm, TYPE-C Cartridges. Hichrom Brand. 1 EA.



Additional Info

For Cogent TYPE-C Brand Semi Prep Guard Columns (Cartridges) See Respective Column / Stationary Phases pages on our website.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$832.05

Printed On 2024/05/06