

## Catalog Number

82000-SP-GH

**Guard Column Holder, Semi Prep, 10mm TYPE-C Guard. Accepts 10mm ID x 10mm, TYPE-C Cartridges. Hichrom Brand. 1 EA.**



### Additional Info

For Cogent TYPE-C Brand Semi Prep Guard Columns (Cartridges) See Respective Column / Stationary Phases pages on our website.

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

**Price**  
**\$832.05**

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$832.05**