

## Catalog Number

82060-P30-250

**Prep HPLC Column, NP Diol, 10um, 100A. 30mm ID x 250mm Length.**  
Cogent Brand. 1 EA.



### Additional Info

This Cogent NP Diol™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. Non end-capped, it is manufactured with the highest purity Silica available and bonded with a short hydrocarbon ligand that has a terminal, Diol Functional group.

The 10µm Particle in this Column works well with higher flow rates and yet produces lower back pressures than smaller particles. The large diameter this column provides good loading capacity which is great for purification of compounds making it ideal for Preparative HPLC.

*Use for hydrophilic or polar compounds in HILIC or Normal Phase Prep HPLC instead of unmodified Silica Columns.*

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$9,152.44**

Printed On 2024/04/29