

Catalog Number

82061-05P-4-3M

HPLC Column, NP Silica, Wide Pore, 10um, 4.0mm ID x 50mm Length, 300A. Cogent Brand. 1 EA. CUSTOM ORDER.



Additional Info

This Cogent NP Silica™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. This low metal content, 300A, very high purity Silica is not modified in any way and is ideal for new methods involving large molecules.

With a 10µm Particle, this column works well with higher flow rates and lower back pressures compared with smaller particles,

Often used as a scouting column for large hydrophilic or polar compounds, in preparative HILIC or Normal Phase HPLC modes.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$576.05

Printed On 2024/05/08