

# Catalog Number

82061-15P-4

**HPLC Column, NP Silica, 10m, 4.0mm ID x 150mm Length, 100A.**  
Cogent Brand. 1 EA. CUSTOM ORDER.



## Additional Info

This Cogent NP Silica™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. This low metal content, 100A, very high purity Silica is not modified in any way and is ideal as an Alternate in any USP Method that uses the L3 Classification or for new methods.

With a 10µm Particle, this column works well with higher flow rates and lower back pressures compared with smaller particles,

*Often used for hydrophilic or polar compounds in HILIC or Normal Phase HPLC.*

**Price**  
**\$674.70**

## Additional Info

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$674.70**