

**Catalog Number****82061-15P-4-3M**

**HPLC Column, NP Silica, Wide Pore**, 10 $\mu$ m, 4.0mm ID x 150mm Length, 300A. Cogent Brand. 1 EA. CUSTOM ORDER.

**Additional Info**

This Cogent NP Silica™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. This low metal content, 300A, very high purity Silica is not modified in any way and is ideal for new methods involving large molecules.

With a 10 $\mu$ m Particle, this column works well with higher flow rates and lower back pressures compared with smaller particles,

*Often used for large hydrophilic or polar compounds in HILIC or Normal Phase HPLC*

**Price**  
**\$674.70**

## Additional Info

*modes.*

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$674.70**