

## Catalog Number

82061-25P

**HPLC Column, NP Silica, 10m, 4.6mm ID x 250mm Length, 100A.**  
Cogent Brand. 1 EA.



### Additional Info

This very Popular Cogent NP Silica™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. This low metal content, 100A, very high purity Silica is not modified in any way and is ideal as an Alternate in any USP Method that uses the L3 Classification or for new methods.

With a 10µm Particle, this column works well with higher flow rates and lower back pressures compared with smaller particles,

*Often used for hydrophilic or polar compounds in HILIC or Normal Phase HPLC.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$733.16**

Printed On 2024/08/06

## Additional Info

**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451  
tel. (732) 380-8900, fax (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$733.16**

Printed On 2024/08/06