

Catalog Number

82061-P21-1003M

Prep HPLC Column, NP Silica, Wide Pore, 10um, 300A. 21.2mm ID x 100mm Length. Cogent Brand. 1 EA.



Additional Info

This Cogent NP Silica™ Wide Pore, Stationary Phase has been a premium value, and used by HPLC laboratories and pilot plants for many years. This low metal content, 100A, very high purity Silica is not modified in any way and is ideal for purification of lipids or other polar compounds.

The 10 µm Particle works well with higher flow rates and produces lower back pressures than smaller particles. The large diameter and length of this column provides high loading capacity which is great for purification of large hydrophilic or charged compounds making this it ideal for Preparative HPLC.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$3,287.15

Printed On 2024/08/06

Additional Info

A specialized HPLC System is required for this columns which is often used in HILIC or Normal Phase mode.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$3,287.15

Printed On 2024/08/06