

Catalog Number

82061-P30-PG1-3

Prep Guard Column Cartridge, NP Silica, Wide Pore, 10um, 300A.
30mm ID x 10mm Length. Cogent Brand. 1 EA.



Additional Info

This Guard Column Cartridge is packed with Wide Pore, 10µm, Cogent NP Silica™ Stationary Phase. This low metal content, 300A, very high purity Silica is not modified in any way.

It is designed to protect your NP Silica Prep Column from particulate and contaminants that will adsorb to it. Packed in a Stainless Steel Tube, it has SS Frits.

Use this Guard Column Cartridge with a 30mm ID Prep Column Holder. Click [HERE](#) for more information.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$537.54

Printed On 2024/08/06

Additional Info

ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$537.54

Printed On 2024/08/06