

Catalog Number

82061-SP250

Semi Prep HPLC Column, NP Silica, 10um, 100A. 10mm ID x 250mm Length. Cogent Brand. 1 EA.



Additional Info

This Cogent NP Silica™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. This low metal content, 100A, very high purity Silica is not modified in any way.

The 10µm Particle works well with higher flow rates and produces lower back pressures than smaller particles. The diameter and length of this column provides increased loading capacity compared with analytical columns, which is why it is chosen for purification of hydrophilic or polar compounds, making it ideal for Semi Preparative HPLC.

A specialized HPLC System is required for this column which is often used in HILIC or Normal

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$2,384.79

Printed On 2024/08/06

Additional Info

Phase mode.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$2,384.79

Printed On 2024/08/06