

## Catalog Number

82361-25P

**HPLC Column, NP Silica, 3um, 4.6mm ID x 250mm Length, 100A.**  
Cogent Brand. 1 EA.



### Additional Info

This Cogent NP Silica™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. This low metal content, 100A, very high purity Silica is not modified in any way and is ideal as an Alternate in any USP Method that uses the L3 Classification or for new methods.

With a 3um Particle, this column is more efficient than a 5um column of the same phase and you will experience sharper peaks. However, due to the smaller particle size, there is an increased back pressure.

*Often used for hydrophilic or polar compounds in HILIC or Normal Phase HPLC.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,036.93**

Printed On 2024/05/08