

Catalog Number

82534-15P-2

HPLC Column, NP Cyano, 5um, 2.1mm ID x 150mm Length, 100A.
Cogent Brand. 1 EA.



Additional Info

This Cogent NP Cyano™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. Not end capped, it is manufactured with a 3 carbon ligand that has an attached terminal Cyano (*Carbon & Nitrogen*) group which includes a triple bond, that is bonded to the highest purity silica available making this Column ideal as an Alternate in USP Methods that use the L10 Classification in Normal Phase.

With a 5 µm Particle, this Column works well with higher flow rates and lower back pressures compared with smaller Particles. It is often used for polar compounds in

Price

\$921.70

Additional Info

Normal Phase HPLC and you should expect shorter retention times and faster equilibration (*between runs*) than with a non-modified, Silica Column.

Packed and shipped in Normal Phase Solvents; click [HERE](#) for Cyano Columns packed in Reversed Phase (polar) Solvents.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$921.70