

## Catalog Number

82561-15P-4-3M

**HPLC Column, NP Silica, Wide Pore, 5µm, 4.0mm ID x 150mm Length, 300A. Cogent Brand. 1 EA. CUSTOM ORDER.**



### Additional Info

This Cogent NP Silica™ Wide Pore Stationary Phase has been a premium value and used by HPLC laboratories for many years. This low metal content, 300A, very high purity Silica is not modified in any way and is ideal as an Alternate in any USP Method that uses the L3 Classification for large molecules or for new methods.

With a 5 µm Particle and 300A Pore Size, this column works well with higher flow rates and lower back pressures.

*Often used for large hydrophilic or polar compounds in HILIC or Normal Phase HPLC.*

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$848.36**

Printed On 2024/04/28