

Catalog Number**82561-25P-4-3M****HPLC Column, NP Silica, Wide Pore, 5 μ m, 4.0mm ID x 250mm Length, 300A. Cogent Brand. 1 EA. CUSTOM ORDER.****Additional Info**

This Cogent NP Silica™ Wide Pore Stationary Phase has been a premium value and used by HPLC laboratories for many years. This low metal content, 300A, very high purity Silica is not modified in any way and is ideal as an Alternate in any USP Method that uses the L3 Classification for large molecules or for new methods.

With a 5 μ m Particle and 300A Pore Size, this column works well with higher flow rates and lower back pressures.

Often used for large hydrophilic or polar compounds in HILIC or Normal Phase HPLC.

Price
\$942.44

Additional Info

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$942.44