

# Catalog Number

82561-PG2-3

**Guard Column Cartridges, NP Silica Wide Pore, 5µm, 300A. 4.0mm ID x 10mm Length. Cogent Brand. 5/PK.**



## Additional Info

These Guard Column Cartridges are packed with Wide Pore, 5µm, Cogent NP Silica™ Stationary Phase. This low metal content, 300A, very high purity Silica is not modified in any way.

Designed to protect your 3.0, 4.0 or 4.6mm ID, NP Silica™, Analytical Columns from particulate and contaminants that will adsorb to it. They are packed in Stainless Steel Tubes and have Stainless Steel Frits that are pressed into PEEK Frit Caps.

*These Guard Column Cartridges must be used with a Guard Column Holder,*

**Price****\$458.03**

## Additional Info

click [HERE](#) for more information about the Holder.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$458.03**