

## Catalog Number

82561-SP250-3M

**Semi Prep HPLC Column, NP Silica Wide Pore, 5um, 300A. 10mm ID x 250mm Length. Cogent Brand. 1 EA.**



### Additional Info

This Cogent NP Silica™ Wide Pore, Stationary Phase has been a premium value, and used by HPLC laboratories and pilot plants for many years. This low metal content, 100A, very high purity Silica is not modified in any way and is ideal for purification of lipids or other polar compounds.

The 5µm Particle works well with prep flow rates and produces lower back pressures than smaller particles. The large diameter and length of this column provides high loading capacity which is great for purification of large hydrophilic or polar compounds making this ideal for Preparative HPLC.

*A specialized HPLC System is required for this columns which is often used in HILIC or Normal Phase mode.*

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$2,674.84**

Printed On 2024/05/05