

Catalog Number

83000-GH

Guard Column Holder, Accepts 1.0-4.0 mm ID, 10 mm long Cogent brand guard cartridges. Cogent Brand. 1 EA.



Additional Info

For Cogent RP, Cogent NP or Cogent IEX Guard Columns (Cartridges) See Respective Column / Stationary Phases pages on this website.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$313.59

Printed On 2024/07/20