

## Catalog Number

## 95001-03ND

**Inserts, Low Volume** for Vials, Glass. Clear, Deactivated, 100ul fill volume with a conical, precision point and 5x30mm outer dimensions with plastic spring attached. For use in 8mm, Push Top, 1ml autosampler



### Additional Info

First hydrolytic, Class A, Type 33, Deactivated / Silanized Borosilicate Glass. Features thicker glass walls than other similar inserts and a precision point conical tip for maximum recovery of sample. Residual volume is 2ul with most autosamplers.

Polypropylene springs are attached to keep the inserts center aligned and reduce chances of the needle breaking the inserts during use. This product is manufactured in an ISO 9001:2015 Certified facility.

Can be used with 8x40mm, "96 position" vials.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$140.63**

Printed On 2024/05/18