

## Catalog Number

95001-0AP

**Vial & Cap Kit** includes 250 EA 1ml, Push Top, Amber, 8x40mm, Autosampler Vials and 250 EA matching Push Top Caps. AQ Brand. 250/PK.



### Additional Info

1st hydrolytic, Type N51A Borosilicate Glass, amberized and considered low actinic. This product is manufactured in an ISO 9001:2015 Certified facility.

Push caps are made with ultra high purity polyethylene and have a thinned penetration point in a star configuration in the cap center to allow venting during injection and prevents the cap from grabbing the needle when it is extracted from the vials.

*Use with Waters® WISP® 96 Position / Carousel autosamplers.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$70.95**

Printed On 2024/08/06

## Additional Info

**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451  
tel. (732) 380-8900, fax (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$70.95**

Printed On 2024/08/06