

## Catalog Number

## 95001-INS

**Inserts, Low Volume** for Vials, Glass. Clear, 150ul fill volume with a flat precision point and 5x34mm outer dimensions. For use in 8mm, Push Top, 1ml autosampler vials. MicroSolv Brand. 1,000/CS.



### Additional Info

First hydrolytic, Class A, Type N51A, Borosilicate Glass. Features a flat precision point conical tip for maximum recovery of sample. Residual volume is 4ul with most autosamplers.

Spring may be required, sold separately, recommended to *use with 95001-OCP-INS caps, in the 95001-OCP 8x40mm, 96 Position vials.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$672.05**

Printed On 2024/05/18