

Catalog Number

95001-INS

Inserts, Low Volume for Vials, Glass. Clear, 150ul fill volume with a flat precision point and 5x34mm outer dimensions. For use in 8mm, Push Top, 1ml autosampler vials. MicroSolv Brand. 1,000/CS.



Additional Info

First hydrolytic, Class A, Type N51A, Borosilicate Glass. Features a flat precision point conical tip for maximum recovery of sample. Residual volume is 4ul with most autosamplers.

Spring may be required, sold separately, recommended to *use with 95001-0CP-INS caps, in the 95001-0CP 8x40mm, 96 Position vials.*

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$672.05

Printed On 2024/07/27