

## Catalog Number

95010-E-0AV-K

**Vials, Headspace, Screw Top, Glass. Amber, 10ml, 18mm precision threads, 22.5x46mm outer dimensions, Rounded Bottom. CTC Compliant. BASIK Brand. (10 packs of 100) 1000/CS.**



### Additional Info

First hydrolytic, Class A, Borosilicate Glass, amberized and considered low actinic. Features 18mm precision threads for proper sealing.

Compatible with CTC®, LEAP Combi Pal®, Agilent® and other Headspace autosamplers. Not suitable for older versions of TurboMatrix 110® that were manufactured before 2006.

Lowest cost in class. Click [HERE](#) for smaller pack quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$887.40**

Printed On 2024/05/18