

Catalog Number

95010-E-0CV

Vials, Headspace, Screw Top, Glass. Clear, 10ml, 18mm precision threads, 22.5x46mm outer dimensions, Rounded Bottom. CTC Compliant. BASIK Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Borosilicate Glass. Features 18mm precision threads for proper sealing. Compatible with CTC®, LEAP Combi Pal®, Agilent® and other Headspace autosamplers. Not suitable for older versions of TurboMatrix 110® that were manufactured before 2006.

Lowest cost in class. Click [HERE](#) for case quantities.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$73.96

Printed On 2024/08/06

Additional Info

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$73.96

Printed On 2024/08/06