

## Catalog Number

95020-21S-T

**Caps, Headspace**, Screw Top, with White Silicone Rubber / Red PTFE Pre-Slit (Star Shape) Septa, fitted in Steel, Magnetic, Silver Caps with 8mm center holes. For use with 18mm Precision Thread Headspace



### Additional Info

Star Shape Pre-Slit provides a greater target area for needles and are easier to penetrate but may not reseal as well as single or non-slit septa. Recommended for SPME autosamplers including CTC Combi Pal®.

Silicone Rubber / UltraClean PTFE Septa is 55° Shore A (hardness) and is 1.5mm thick. Cap has ND18 Screw Threads.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$224.15**

Printed On 2024/05/07