

## Catalog Number

95020-26-S

**Caps, Headspace, Screw Top, with White Silicone Rubber / Blue PTFE SPME Style Septa, fitted in Steel, Magnetic, Silver Caps with 8mm center holes. For use with 18mm Precision Thread Headspace Vials.**



### Additional Info

Special Shape Septa with thick outer ring and thinned center provides a good sealing with the vial and a target area easier for needles to penetrate. Recommended with CTC Combi Pal® and other autosamplers with SPME or Solid Phase Micro-Extraction injectors.

Cap has ND18 Screw Threads.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$84.60**

Printed On 2024/08/06

## Additional Info

**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451  
tel. (732) 380-8900, fax (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$84.60**

Printed On 2024/08/06