

## Catalog Number

95020-50S-W

**Caps, Headspace**, Screw Top, with White Silicone Rubber / Blue PTFE Septa, fitted in Steel, Magnetic, Silver Caps with 8mm center holes. For use with 18mm Precision Thread Headspace Vials. BASIK Brand.



### Additional Info

Silicone Rubber / PTFE Septa has durometer of 50° Shore A (hardness) and is 1.5mm thick. Cap has ND18 Screw Threads. Click [HERE](#) for case quantities.

Economical and low cost for budget minded laboratories.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$120.31**

Printed On 2024/05/06