

Catalog Number

95020-E-0CV-K

Vials, Headspace, Screw Top, Glass. Clear, 20ml, 18mm precision threads, 22.5x75.5mm outer dimensions, Rounded Bottom. CTC Compliant. BASIK Brand. (10 packs of 100) 100/PK.



Additional Info

First hydrolytic, Class A, Borosilicate Glass, Features 18mm precision threads for proper sealing. Compatible with CTC®, LEAP Combi Pal®, Agilent® and other Headspace autosamplers. Not suitable for older versions of TurboMatrix 110® that were manufactured before 2006.

Lowest cost in class. Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$845.13

Printed On 2024/04/29