

Catalog Number

95025-01-3D

Caps, Headspace, Crimp Top, with Dark Grey Butyl Rubber Septa, fitted in Steel, Magnetic, Gold Caps with 5mm center holes. For use with 20mm Crimp Top Headspace Vials. MicroSolv Brand. 100/PK.



Additional Info

Molded, Butyl Rubber Septa is 55° Shore A (hardness) and is 3.0mm thick. Temperature resistant from -40C to 120C. Recommended for non critical analyses or when PTFE is an analyte or must be avoided. Recommended for CTC™ and other autosamplers.

Low cost for budget minded laboratories.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$50.40

Printed On 2024/05/18