


## Catalog Number

95025-02-2S

**Caps, Headspace, Crimp Top, with Grey Butyl Rubber Laminated PTFE Septa, fitted in Aluminum, Silver Pressure Relief Caps with 8mm center holes. For use with 20mm Crimp Top Headspace Vials. MicroSolv Brand.**



### Additional Info

Butyl Rubber Septa that is laminated with PTFE on both sides and is 55° Shore A (hardness) and is 3.0mm thick. Temperature resistant from -40C to 120C.

Clear lacquered Caps have score lines that will break open at 3.0 +/- 0.5 Bar, for safety. Low cost for budget minded laboratories.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$94.31**

Printed On 2024/08/06

## Additional Info

**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451  
tel. (732) 380-8900, fax (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$94.31**

Printed On 2024/08/06