


## Catalog Number

95025-04

**Septa, 20mm, Headspace, Blue Silicone Rubber completely laminated with Transparent PTFE.** For use in 20mm Crimp Caps, sold separately. MicroSolv Brand. 1000/CS.



### Additional Info

Silicone Rubber with Durometer of 45 Shore A, (hardness) and 3.0mm thick. Completely laminated with Ultra Clean PTFE on both sides of the Septa which is good for sensitive analyses. The lower Durometer means that this is a softer Septa and therefore better for easier penetrations by autosampler needles. Wide temperature range from -60° to 200°C.

Use in Aluminum or Steel Crimp Caps for Headspace autosamplers. Supplied in Polyethylene Bags.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,206.94**

Printed On 2024/05/07