


Catalog Number

95025-04-1S

Caps, Headspace, Crimp Top, with Translucent Blue Silicone Rubber / White PTFE Septa, fitted in Aluminum, Silver Caps with 10mm center holes. For use with 20mm Crimp Top Headspace Vials. MicroSolv Brand.



Additional Info

Silicone Rubber / Ultraclean PTFE with Durometer of 45 Shore A, (hardness) and 3.0mm thick. Completely laminated with PTFE on both sides of the Septa which is good for sensitive analyses. The lower Durometer means that this is a softer Septa and therefore better for easier penetrations by autosampler needles. Wide temperature range from -60° to 200°C.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$138.71

Printed On 2024/05/01