

Catalog Number

95025-04-5R

Caps, Headspace, Crimp Top, with Transparent Blue Silicone Rubber / Transparent PTFE Septa, fitted in Bimetallic, Magnetic, Red Caps with 8mm center holes. For use with 20mm Crimp Top Headspace Vials.



Additional Info

Silicone Rubber / Ultraclean PTFE with Durometer of 45 Shore A, (hardness) and 3.2mm thick. Completely laminated with PTFE on both sides of the Septa which is good for sensitive analyses. The lower Durometer means that this is a softer Septa and therefore better for easier penetrations by autosampler needles. Wide temperature range from -60° to 200°C.

Bimetallic caps are still magnetic but easier to crimp compared to all steel crimp caps.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$287.11

Printed On 2024/05/04