

Catalog Number

95025-08-5B

Caps, Headspace, Crimp Top, with White Silicone Rubber / White PTFE, High Purity Septa, fitted in Bimetallic, Magnetic, Blue Caps with 8mm center holes. For use with 20mm Crimp Top Headspace Vials. BASIK



Additional Info

SIMILAR TO 95025-09-2S

High Purity Silicone Rubber / PTFE Septa with Durometer of 45 Shore A, (hardness) and 3.0mm thick with Ultra Low Bleed Rubber that is made with a proprietary process for ultra durability, purity and low extractables.

Excellent choice for highly sensitive analyses. The lower Durometer means that this is a softer Septa and therefore better for easier penetrations by autosampler needles.

Wide temperature range from -60_x0081_° to 200_x0081_°C.

Price
\$83.00

Additional Info

Bimetallic caps are still magnetic but easier to crimp compared to all steel crimp caps.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$83.00