


Catalog Number

95025-09

Septa, 20mm, Headspace, Translucent Silicone Rubber / Natural PTFE.
For use in 20mm Crimp Caps, sold separately. Ultra Low Bleed. BASIK Brand. 1000/CS.



Additional Info

Silicone Rubber / PTFE Septa with Durometer of 45 Shore A, (hardness) and 3.0mm thick with Ultra Clean PTFE and Ultra Low Bleed Rubber that is made by a proprietary process for ultra high purity, low extractables and temperature durability.

This Septa is good for highly sensitive analyses. No color has been added to this Septa. The lower Durometer means that this is a softer Septa and therefore better for easier penetrations by autosampler needles. Wide temperature range from -60° to 200°C.

Use in Aluminum or Steel Crimp Caps for Headspace autosamplers. Supplied in Polyethylene Bags.

This is a Custom Order Item. Please contact Customer Service for Details.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$499.01

Printed On 2024/05/04