

Catalog Number

95025-11-1S

Caps, Headspace, Crimp Top, with Natural Silicone Rubber / Natural PTFE (10mil), High Purity Septa, fitted in Aluminum, Silver Caps with 10mm center holes. For use with 20mm Crimp Top Headspace Vials.



Additional Info

High Purity Silicone Rubber / PTFE Septa (10mil thick) with Durometer of 50 Shore A, (hardness) and 2.5mm thick.

Wide temperature range from -60° to 200°C.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$53.70

Printed On 2024/05/05