

Catalog Number

95025-11-2-1S

Caps, Headspace, Crimp Top, with Natural Silicone Rubber / Natural PTFE (5mil), High Purity Septa, fitted in Aluminum, Silver Caps with 10mm center holes. For use with 20mm Crimp Top Headspace Vials.



Additional Info

High Purity Silicone Rubber / PTFE Septa (5mil thick) with Durometer of 50 Shore A, (hardness) and 3.0mm thick.

Wide temperature range from -60° to 200°C.

Price
\$53.52

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$53.52