

Catalog Number

95025-13-1S

Caps, Headspace, Crimp Top, with White Silicone Rubber / Ivory PTFE (3mil), High Purity Septa, fitted in Aluminum, Silver Caps with 10mm center holes. For use with 20mm Crimp Top Headspace Vials. BASIK



Additional Info

High Purity Silicone Rubber / PTFE Septa (3mil thick) with Durometer of 45 Shore A, (hardness) and 3.0mm thick.

Wide temperature range from -60° to 200°C.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$53.68

Printed On 2024/05/04