

## Catalog Number

95025-14

**Septa, 20mm, Headspace, White Silicone Rubber / Silver Foil.** For use in 20mm Crimp Caps, sold separately. High Temperature. MicroSolv Brand. (10x100) 1000/CS.



### Additional Info

Silicone Rubber / Aluminum Foil Silver Septa with Durometer of 50 Shore A, (hardness) and 3.0mm thick is good for high temperature analyses. No color has been added to this Septa. Often used with Perkin Elmer Instruments.

Wide temperature range from -60° to 220°C.

Use in Aluminum or Steel Crimp Caps for Headspace autosamplers. Supplied in Polyethylene Bags.

**Price**

**\$2,049.71**

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$2,049.71**