



# Catalog Number

# 95025-15-1S

**Caps, Headspace**, Crimp Top, with Orange Silicone Rubber / Foil, fitted in Aluminum, Silver Caps with 10mm center holes. For use with 20mm Crimp Top Headspace Vials. BASIK Brand. 100/PK.



## Additional Info

Silicone Rubber / Aluminum Foil Silver Septa with Durometer of 50 Shore A, (hardness) and 2.5mm thick is good for high temperature analyses. No color has been added to this Septa. Often used with Perkin Elmer Instruments.

Wide temperature range from -60° to 220°C.

**Price**  
**\$71.21**



## Additional Info

### **MICROSOLV**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### **Price**

**\$71.21**