


Catalog Number

95025-23-4D

Caps, Headspace, Crimp Top, with White Silicone Rubber / Blue PTFE Septa, fitted in Steel, Magnetic, Gold Caps with 8mm center holes. For use with 20mm Crimp Top Headspace Vials. MicroSolv Brand. 100/PK.



Additional Info

Compatible with SPME, Solid Phase Micro-Extraction capable autosamplers including CTC Pal®.

Silicone Rubber / UltraClean PTFE Septa is 55° Shore A (hardness) and is 1.5mm thick. Use with SPME Vials. Click [HERE](#) to view SPME Vials.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$181.74

Printed On 2024/04/29