

Catalog Number

95025-26-4B

Caps, Headspace, Crimp Top, with White Silicone Rubber / Blue PTFE SPME Style Septa, fitted in Bimetallic, Magnetic, Blue Caps with 8mm center holes. For use with 20mm Crimp Top Headspace Vials. BASIK



Additional Info

Special Shape Septa with thick rubber outer ring (3mm) and thinned rubber center (1.5mm) provides a good seal with most 20mm crimp top Headspace vials. The target area on this Septa is easier for SPME needles to penetrate, minimizing broken or damaged needles.

Recommended with CTC Pal® and other autosamplers with SPME or Solid Phase Micro-Extraction injectors.

Bimetallic Caps have a magnetic center and an aluminum outer ring that makes crimping them easier than all steel caps.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$103.99

Printed On 2024/08/06

Additional Info

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$103.99

Printed On 2024/08/06