


## Catalog Number

95025-60-C

**Caps, Headspace**, Crimp Top, with Grey Butyl Rubber / Transparent PTFE Septa, fitted in Steel, Magnetic, Gold Caps with 8mm center holes. For use with 20mm Crimp Top Headspace Vials. BASIK Brand. 1000/CS.



### Additional Info

Butyl Rubber Septa that is laminated with PTFE and is 55<sup>0</sup> Shore A (hardness) and is 3.0mm thick.

Temperature resistant from -40C to 120C.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,494.83**

Printed On 2024/05/03