

Catalog Number

95025-AG-06C

Vials, Headspace, Crimp Top, Glass. Clear, 6ml, 20mm Serum Finish, 22x38mm outer dimensions, Flat Bottom with a Long Neck. Agilent Compatible. BASIK Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Borosilicate Glass, Use with aluminum or steel crimp caps for proper sealing. Compatible with Agilent®, Carlo Erba®, Dani® and other Headspace autosamplers.

These vials are also used with Karl Fisher titrators by Mettler and Metrohm.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$63.28

Printed On 2024/04/28