

Catalog Number

95025-AG-06C-D

Vials, Headspace, Crimp Top, Glass. Clear, 6ml, Deactivated, 20mm Serum Finish, 22x38mm outer dimensions, Flat Bottom with a Long Neck. Agilent Compatible. BASIK Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Deactivated / Silanized Borosilicate Glass, Use with aluminum or steel crimp caps for proper sealing. Compatible with Agilent®, Carlo Erba®, Dani® and other Headspace autosamplers.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$133.20

Printed On 2024/04/29