


Catalog Number

95025-E-AG10A-K

Vials, Headspace, Crimp Top, Glass. Amber, 10ml, 20mm Serum Finish, 22.5x46mm outer dimensions, Flat Bottom with a Long Neck. Agilent Compatible. BASIK Brand. 1000/CS.



Additional Info

First hydrolytic, Class A, Borosilicate Glass, amberized and considered low actinic. Use with aluminum or steel crimp caps for proper sealing.

Compatible with Agilent® and other Headspace autosamplers. Click [HERE](#) for smaller pack quantities.

Price
\$728.25

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$728.25