

Catalog Number

95025-E-AG20C

Vials, Headspace, Crimp Top, Glass. Clear, 20ml, 20mm Serum Finish, 23x75.5mm outer dimensions, Flat Bottom with a Long Neck. Agilent Compatible. BASIK Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Borosilicate Glass, Use with aluminum or steel crimp caps for proper sealing. Compatible with Agilent®, Carlo Erba®, Dani® and other Headspace autosamplers.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$71.64

Printed On 2024/05/02