


Catalog Number

95025-E-AP06A

Vials, Headspace, Crimp Top, Glass. Amber. 6ml, 20mm Serum Finish, 22x38mm outer dimensions, Beveled Top and Flat Bottom with a Tapered Neck. Agilent Compatible. BASIK Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Borosilicate Glass, amberized and considered low actinic. Use with aluminum or steel crimp caps for proper sealing.

Compatible with Agilent®, Perkin Elmer® and other Headspace autosamplers.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$66.95

Printed On 2024/04/27